

Minimum Volume IV Extension Sets

Our Minimum Volume IV Extension Sets have an inner diameter of 1.5mm and are available in a range of lengths.

- ⊗ Transparent, soft, smooth tubing
- ⊗ 1.5mm inner diameter
- ⊗ Female Luer lock to male Luer lock with rotating collar
- ⊗ Made of non-toxic medical grade PVC



Minimum Volume IV Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV051002	Minimum Volume Extension Set with Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	600cm	—	50
IV051003	Minimum Volume Extension Set with Injection Site Female Luer Lock to Male Luer Lock & RC	165cm	Standard Injection Site	100

Table continued on next page

Minimum Volume IV Extension Sets

Minimum Volume IV Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV051004	Minimum Volume Extension Set with Needleless Access Site Female Luer Lock to Male Luer Lock and RC	165cm	Needleless Access Site	100
IV051006	Minimum Volume Extension Set with Needleless Access Site, Non-Return valve, Female Luer Lock to Male Luer Lock and RC	20cm	Needleless Access Site	100
IV051010	Minimum Volume Extension Set with Anti-Siphon Valve	150cm	—	100
IV054000	Minimum Volume Extension Set with 1 single lumen with Needleless Access Site and Single Lumen Clamp	15cm	Needleless Access Site	200