


Standard Bore Extension Sets

Our Standard Bore Extension Sets have an inner diameter of 3mm and are available in a range of lengths and configurations.

- ⊙ Transparent, soft, smooth tubing
- ⊙ 3mm inner diameter
- ⊙ Female Luer lock to male Luer lock with rotating collar
- ⊙ Made of non-toxic medical grade PVC



Standard Bore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV052000	Standard Bore Extension Set with Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	75cm	—	100
IV052001	Standard Bore Extension Set with Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	150cm	—	100

Table continued on next page

Standard Bore Extension Sets

Standard Bore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV052100	Standard Bore Extension Set - Injection Site Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	495cm	Standard Injection Site	50
IV052101	Standard Bore Extension Set – Needleless access site, Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	495cm	Needleless Access Site	50
IV052201	Standard Bore Extension Set - 2 x Needleless Access Site Female Luer Lock to Male Luer Lock and RC	80cm	2 x Needleless Access Site	100
IV052203	Standard Bore Extension Set with Female Luer Lock to Male Luer Lock and 2 x Luer Activated Valves	120cm	2 x Luer Activated Valves	100