


Vented Spike Adaptor

Designed for safe and efficient fluid transfer, the m|devices Vented Spike Adaptor enables the transfer of fluids from bags, vials or bottles without the need for needles. The needle-free access site minimises the risk of needle stick injuries, promoting user safety.

- ⊗ Vented spike
- ⊗ Standard bore tubing 100mm in length
- ⊗ Click clamp
- ⊗ Needleless access site



Vented Spike Adaptor				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV070004	Vented Spike Adaptor	10cm	Needleless Access Site	120