

Double Spike Adaptor

Our 30cm Double Spike Adaptor is designed with two vented spikes with a roller regulator and a septum connector.

- ⊗ Allows adaption from single bag to double bag
- ⊗ Incorporates a sealed septum
- ⊗ 30cm length
- ⊗ Free from latex
- ⊗ Made from medical grade, non-toxic PVC



Double Spike Adaptor

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
IV071000	Infusion/Transfusion Double Spike Adaptor	30cm	100