


## Vented Double Spike Adaptor

Our 30cm Vented Double Spike Adaptor is designed with two vented spikes, two roller regulators and a sealed septum connector.

- ⊗ Allows adaption from single bag to double bag
- ⊗ Incorporates a sealed septum
- ⊗ 30cm length
- ⊗ Not made with natural rubber latex
- ⊗ Made from medical grade, non-toxic PVC



### Vented Double Spike Adaptor

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
IV071001	Infusion/Transfusion Vented Double Spike Adaptor	30cm	50