


Closed System Transfer Set

Our Closed System Transfer Set is a needle free device with multi-lumen extensions designed for efficient fluid transfer

It consists of 2 x 2-lumen female luer lock connectors, attached to a 3-way stopcock and 3 perforators with click clamps and a swabable needleless access valve.

- ⊗ 43cm total length
- ⊗ 3mm diameter tubing
- ⊗ 2 x 2-lumen set
- ⊗ 3 spike perforators
- ⊗ Controlled by a 3-way stopcock



Closed System Transfer Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
IV072001	Closed System Transfer Set	43cm	20

m|devices[®]

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321