


Stopcocks (Luer Access)

Our 3-Way stopcocks are available as standard and lipid resistant versions. They feature transparent channels and one rotating male Luer connector as well as two female Luer connectors. Arrows on the handle indicate the direction of flow.

- ⊗ 3-Way stopcock
- ⊗ For pressure up to 65psi
- ⊗ Male Luer lock with rotating collar
- ⊗ Two female Luer connections
- ⊗ Standard and lipid resistant versions
- ⊗ Flow indication arrows



| 3-Way Stopcocks with Luer Access | | | STERILE |
|----------------------------------|----------------------------------|--------|------------|
| RE-ORDER | DESCRIPTION | COLOUR | CARTON QTY |
| IV110000 | Stopcock - 3-Way | White | 50 |
| IV110001 | Stopcock - 3-Way Lipid Resistant | Blue | 50 |
| IV110011 | Stopcock - 3-Way Lipid Resistant | Yellow | 50 |