

## Stopcocks (Extension Tubing)

Our Stopcock Extension Sets are made from medical grade, non-toxic, kink resistant 3mm diameter tubing and are available in a range of lengths.

- ⊗ 3-way stopcock incorporating extension tubing
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar
- ⊗ Two female Luer connectors
- ⊗ Flow indicator arrows



### 3-Way Stopcocks with Extension Tubing

STERILE

RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
IV110002	Stopcock Lipid Resistant with 15cm Extension Set	15cm	Blue	50
IV110003	Stopcock Lipid Resistant with 25cm Extension Set	25cm	Blue	50
IV110004*	Stopcock Lipid Resistant with 50cm Extension Set	50cm	Blue	50
IV110005*	Stopcock Lipid Resistant with 80cm Extension Set	80cm	Blue	50
IV110006*	Stopcock Lipid Resistant with 100cm Extension Set	100cm	Blue	50

\*Minimum order quantity required