

3-Way Stopcock

Our lipid resistant 3-Way Stopcock is simple to use, featuring a handle that easily rotates into position. It also requires a minimum amount of priming fluid.

- ⊗ Lipid resistant
- ⊗ Male Luer lock with rotating collar
- ⊗ Two female Luer connections
- ⊗ Rated up to 450kpa water pressure or 300kpa air pressure



3-Way Stopcocks

STERILE

RE-ORDER	DESCRIPTION	COLOUR	CARTON QTY
IV110011	Stopcock - 3-Way Lipid Resistant	Yellow	50