

## Insulin Pen Needle

Our Insulin Pen Needles consist of: a needle holder, needle tube, tip sheath and needle seat sheath. Available in a range of gauges and lengths

This product can be used with the following brands of insulin pens: Novopen, Wanbangpen, Prandilin, HumaPen, Gansulinpen and KlikSTAR.

- ⊗ Colour coded for easy size identification
- ⊗ Gauge range from 29G-32G
- ⊗ Free from Latex and DEP
- ⊗ Lengths range from 4mm -12.7mm
- ⊗ Sterile EO
- ⊗ Compatible with various popular brands of insulin pens



Insulin Pen Needles				STERILE
RE-ORDER	DESCRIPTION	COLOUR	SIZE	CARTON QTY
IV140531	Insulin Pen Needle 5mm	Purple	31G	100
IV141030	Insulin Pen Needle 10mm	Lime Green	30G	100