

Tuberculin Syringe with Retractable Safety Needle

Our Tuberculin Syringe with Retractable Safety Needle retracts the needle into the barrel after injection to prevent re-use and needle stick injury.

- ⊗ Permanently disabled needle reduces risk of needle stick injury
- ⊗ Eliminates blood splatter
- ⊗ Audible lock
- ⊗ Snap off plunger
- ⊗ Sharps waste reduction



Retractable Safety Syringe

STERILE

RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
IV160113	Tuberculin Syringe with Retractable Safety Fixed Needle 27G x 13mm	1mL	100