

Luer Slip Syringes

Our Luer Slip three-piece Hypodermic syringes are made from Polypropylene, are free from latex and DEHP and are available in a broad range of sizes.

- ⊗ Plunger stopper with double sealing ring
- ⊗ Lipid resistant
- ⊗ Free from latex and DEHP



⊗
Nozzle Positions



Concentric Nozzle

Luer Slip Syringes				STERILE
RE-ORDER	DESCRIPTION	VOLUME	NOZZLE POSITION	CARTON QTY
IV180001	Luer Slip Syringe without Needle	1mL	Concentric	100
IV180003	Luer Slip Syringe without Needle	3mL	Concentric	100
IV180005	Luer Slip Syringe without Needle	5mL	Concentric	100

Table continued on next page