

Luer Slip Syringes - Eccentric

Our Luer Slip three-piece Hypodermic syringes are made from Polypropylene, are free from latex and DEHP. Our eccentric range is available from 20mL upwards.

- ⊙ Flat plunger stopper with double sealing ring
- ⊙ Lipid resistant
- ⊙ Free from latex and DEHP



Nozzle Position



Eccentric Nozzle

Luer Slip Syringes				STERILE
RE-ORDER	DESCRIPTION	VOLUME	NOZZLE POSITION	CARTON QTY
IV180020	Luer Slip Syringe without Needle	20mL	Eccentric	100
IV180030	Luer Slip Syringe without Needle	30mL	Eccentric	75
IV180050	Luer Slip Syringe without Needle	50mL	Eccentric	40
IV180060	Luer Slip Syringe without Needle	60mL	Eccentric	40