

## Luer Slip Syringes - Eccentric

Our Luer Slip three-piece Hypodermic syringes are made from Polypropylene, are free from latex and DEHP. Our eccentric range is available from 20mL upwards.

- ⊗ Flat plunger stopper with double sealing ring
- ⊗ Lipid resistant
- ⊗ Not made with Latex and DEHP



### Nozzle Position



### Eccentric Nozzle

Luer Slip Syringes - Eccentric				STERILE
RE-ORDER	DESCRIPTION	VOLUME	NOZZLE POSITION	CARTON QTY
IV180020	Luer Slip Syringe without Needle	20mL	Eccentric	100
IV180050	Luer Slip Syringe without Needle	50mL	Eccentric	40