

Hypodermic Syringes (Luer Lock)

Our Luer Lock three-piece Hypodermic syringes are made from Polypropylene, are free from latex and DEHP and are available in a broad range of sizes. A retractable safety syringe version retracts the needle into the barrel after injection to prevent re-use and needle stick injury.



- ⊗ Plunger stopper with double sealing ring
- ⊗ Lipid resistant
- ⊗ Free from latex and DEHP

The retractable syringe:

- ⊗ Permanently disabled needle reduces risk of needle stick injury
- ⊗ Eliminates blood splatter
- ⊗ Audible lock
- ⊗ Snap off plunger
- ⊗ Sharps waste reduction

Luer Lock Retractable Safety Syringes			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
IV170003	Luer Lock Syringe - Retractable without Needle	3mL	100
IV170005	Luer Lock Syringe - Retractable without Needle	5mL	100
IV170010	Luer Lock Syringe - Retractable without Needle	10mL	100
IV170020	Luer Lock Syringe - Retractable without Needle	20mL	50

Table continued on next page

Hypodermic Syringes (Luer Lock)

Hypodermic Syringes (Luer Lock)			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
IV181001	Luer Lock Syringe without Needle	1mL	100
IV181003	Luer Lock Syringe without Needle	3mL	100
IV181005	Luer Lock Syringe without Needle	5mL	100
IV181010	Luer Lock Syringe without Needle	10mL	100
IV181020	Luer Lock Syringe without Needle	20mL	100
IV181030	Luer Lock Syringe without Needle	30mL	75
IV181050	Luer Lock Syringe without Needle	50mL	40
IV181060	Luer Lock Syringe without Needle	60mL	40