


Catheter Tip Syringe

Our Catheter Tip Syringe has a specifically designed funnel end that can be connected to a Foley Catheter, to irrigate the catheter as clinically indicated.

- ⊗ Ideal for irrigating catheters and other tubes
- ⊗ Clear bold incremental markings
- ⊗ 60mL capacity
- ⊗ Plunger stopper with double sealing ring
- ⊗ Lipid resistant
- ⊗ Free from latex and DEHP



Catheter Tip Syringe			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
IV182060	Catheter Tip Syringe	60mL	50