

## Luer Activated Valve - Needleless Access Site

Safsite® features a closed silicone diaphragm that checks flow until mechanically opened with a luer device. It allows for bi-directional flow up to 400mL/min upon actuation.

Safsite® is designed for use with luer lock and luer slip connectors. One small push ensures a firm connection also with luer slip syringes. Its design ensures a low priming volume of 0.12mL and a pressure rating of 45 Psi.

- ⊗ Easy multiple access with either luer slip or luer lock syringes
- ⊗ No twist necessary for luer slip connection
- ⊗ Low priming volume
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent housing for clear visibility



Luer Activated Valve			STERILE
RE-ORDER	DESCRIPTION	COLOUR	CARTON QTY
IV200000	Luer Activated Valve	Blue	50