

Luer Activated Valve - Needleless Access Site

This valve has a closed silicone diaphragm that checks flow until mechanically opened with a luer device. It allows for bi-directional flow of up to 400 mL/min upon actuation.

Safsite® is designed for use with luer lock and luer slip connectors. One small push ensures a firm connection also with luer slip syringes. Its design ensures a low priming volume of 0.12mL and a pressure rating of 45 Psi.

- ⊗ Easy multiple access with either luer slip or luer lock syringes
- ⊗ No twist necessary for luer slip connection
- ⊗ Low priming volume
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent housing for clear visibility



Luer Activated Valve

STERILE

RE-ORDER	DESCRIPTION	COLOUR	CARTON QTY
IV200000	Luer Activated Valve	Blue	50

m|devices®

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321