

LuerSite Valve

The LuerSite is designed for use with Luer lock and Luer slip connectors. One small push ensures a firm connection even with Luer slip syringes. A closed silicone diaphragm stops flow until mechanically opened with a Luer device. LuerSite allows for bi-directional flow of up to 345mL/min upon actuation.

It is designed to ensure a low priming volume of 0.14mL and a pressure rating of 45Psi.

- ⊗ Easy, multiple access with either Luer lock or Luer slip connectors
- ⊗ Low priming volume
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent for clear visibility



LuerSite Valve

STERILE

RE-ORDER	DESCRIPTION	COLOUR	CARTON QTY
IV200001	LuerSite Valve - Needleless Access Site	Green	50

m|devices[®]

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321