

LuerSite Needleless Access Site

The LuerSite is designed for use with luer lock and luer slip connectors. One small push ensures a firm connection even with luer slip syringes. A closed silicone diaphragm stops flow until mechanically opened with a luer device. LuerSite allows for bi-directional flow of up to 345mL/min upon actuation.

It is designed to ensure a low priming volume of 0.14mL and a pressure rating of 45 Psi.

- ⊗ Easy multiple access with either luer slip or luer lock syringes
- ⊗ Low priming volume
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent for clear visibility



LuerSite Needleless Access Site			STERILE
RE-ORDER	DESCRIPTION	COLOUR	CARTON QTY
IV200001	LuerSite Valve - Needleless Access Site	Green	50