

Neutral Needleless Access Valve

m|devices Neutral Needleless Access Valve helps to prevent catheter occlusions by reducing blood reflux upon disconnection, in addition to reducing catheter related blood stream infections (CRBSI) through its flat swabable pre-slit septum which allows for easy disinfection.

- ⊗ Flat swabable pre-slit septum
- ⊗ Straight internal fluid pathway
- ⊗ Low priming volume
- ⊗ Gapless space between septum and casing providing a microbial barrier against contamination
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent housing for clear visibility



Neutral Needleless Access Valve

STERILE

RE-ORDER	DESCRIPTION	CARTON QTY
IV200004	Neutral Needleless Access Valve	50

m|devices®

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321