

NeutralSite™

Our NeutralSite™ valve helps to prevent catheter occlusions by reducing blood reflux upon disconnection, in addition to preventing catheter related blood stream infections (CRBSI). Its flat swabable pre-slit septum allows for easy disinfection. In addition, the internal straight fluid pathway helps to minimise the priming volume and address dead space concerns. The NeutralSite™ is recommended for use up to 7 days.

- ⊗ Flat swabable pre-slit septum
- ⊗ Straight internal fluid pathway
- ⊗ Low priming volume
- ⊗ Gapless space between septum and casing providing a microbial barrier against contamination
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent housing for clear visibility
- ⊗ Compatible with luer lock and luer slip connectors



NeutralSite™		STERILE
RE-ORDER	DESCRIPTION	CARTON QTY
IV200005	NeutralSite™	100