

Luer Slip Syringe without Needle


The m|devices Luer Slip Syringes without Needle are intended to draw-up measured medication/solution to be delivered to patients and for the aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.

This range of syringes comes in a variety of sizes from 1mL through to 60mL. They feature a highly transparent barrel with clear, high-contrast graduation, which supports effortless reading and precise dosing.

The conical shaped plunger stopper, with double sealing ring, assists with controlled aspiration and injection and allows easy elimination of air bubbles.

They also features a retaining ring to prevent accidental plunger withdrawal from the barrel.

- ④ Conical shaped plunger stopper to assist with eliminating air bubbles from the syringe
- ④ Highly transparent barrel with high contrast incremental measurement markings
- ④ Double sealing ring for controlled aspiration
- ④ Wide variety of sizes available to meet all clinical requirements



Nozzle Position



Concentric Nozzle

m|devices[®]

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321

Luer Slip Syringe without Needle

Luer Slip Syringe without Needle				STERILE
RE-ORDER	DESCRIPTION	SIZE	NOZZLE POSITION	CARTON QTY
IV480001	Luer Slip Syringe without Needle	1mL	Concentric	100
IV480003	Luer Slip Syringe without Needle	3mL	Concentric	100
IV480005	Luer Slip Syringe without Needle	5mL	Concentric	100
IV480010	Luer Slip Syringe without Needle	10mL	Concentric	100
IV480020	Luer Slip Syringe without Needle	20mL	Concentric	100
IV480030	Luer Slip Syringe without Needle	30mL	Concentric	75
IV480050	Luer Slip Syringe without Needle	50mL	Concentric	40
IV480060	Luer Slip Syringe without Needle	60mL	Concentric	40