


Insulin Pen Needles 3 Bevel

m|devices Insulin Pen Needles 3 Bevel have been designed with a super fine needle to provide a more comfortable injection experience.



The silicone lubricated coating on our triple bevel Insulin Pen Needles provides optimal perforation and easy penetration.

They are available in a range of gauges and lengths to support your clinical requirements.



- ⊗ Sterile packaging
- ⊗ Lubricated needle for a smooth injection
- ⊗ Gauge range: 31-34G
- ⊗ Needle bevel: 3
- ⊗ Needle Lengths available: 4mm, 5mm, 6mm and 8mm
- ⊗ Compatible with various popular brands of insulin pens, refer to IFU for complete listing.

Our Insulin Pen Needles consist of seal, needle hub, needle tube, needle shield and needle container.



Insulin Pen Needles 3 Bevel

Insulin Pen Needles 3 Bevel			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
IV493315	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 5mm	100
IV493316	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 6mm	100
IV493318	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 8mm	100
IV493324	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 4mm	100
IV493325	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 5mm	100
IV493326	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 6mm	100
IV493334	Insulin Pen Needle 3 Bevel	0.2mm (33G) x 4mm	100
IV493344	Insulin Pen Needle 3 Bevel	0.18mm (34G) x 4mm	100