


## Insulin Pen Needles 3 Bevel


m|devices Insulin Pen Needles 3 Bevel have been designed with a super fine needle to provide a more comfortable injection experience.

The silicone lubricated coating on our triple bevel Insulin Pen Needles provides optimal perforation and easy penetration. They are available in a range of gauges and lengths to support your clinical requirements.

- ⊗ Sterile packaging
- ⊗ Lubricated needle for a smooth injection
- ⊗ Gauge range: 31-34G
- ⊗ Needle bevel: 3
- ⊗ Needle Lengths available: 4mm, 5mm, 6mm and 8mm
- ⊗ Compatible with various popular brands of insulin pens, refer to IFU for complete listing.



Our Insulin Pen Needles consist of seal, needle hub, needle tube, needle shield and needle container.



## Insulin Pen Needles 3 Bevel

Insulin Pen Needles 3 Bevel			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
IV493315	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 5mm	100
IV493316	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 6mm	100
IV493318	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 8mm	100
IV493324	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 4mm	100
IV493325	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 5mm	100
IV493326	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 6mm	100
IV493334	Insulin Pen Needle 3 Bevel	0.2mm (33G) x 4mm	100
IV493344	Insulin Pen Needle 3 Bevel	0.18mm (34G) x 4mm	100