


## Insulin Pen Needles 5 Bevel

m|devices, Insulin Pen Needles - 5 Bevel, have been designed with a super fine needle to provide a more comfortable injection experience.


The silicone lubricated coating on our 5 Bevel Insulin Pen Needles provides optimal perforation and easy penetration. Their thin wall technology means less pressing force is required and reduced injection discomfort. It also dramatically reduces the risk of insulin leakage.

They are available in a range of gauges and lengths to support your clinical requirements.

- ⊗ Sterile packaging
- ⊗ Lubricated needle for a smooth injection
- ⊗ Gauge range: 31-34G
- ⊗ Needle bevel: 5
- ⊗ Needle Lengths available: 4mm, 5mm, 6mm and 8mm
- ⊗ Compatible with various popular brands of insulin pens, refer to IFU for complete listing.



Our Insulin Pen Needles consist of seal, needle hub, needle tube, needle shield and needle container.



## Insulin Pen Needles 5 Bevel

Insulin Pen Needles 5 Bevel			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
<b>IV505315</b>	Insulin Pen Needle 5 Bevel	0.25mm (31G) x 5mm	100
<b>IV505316</b>	Insulin Pen Needle 5 Bevel	0.25mm (31G) x 6mm	100
<b>IV505318</b>	Insulin Pen Needle 5 Bevel	0.25mm (31G) x 8mm	100
<b>IV505324</b>	Insulin Pen Needle 5 Bevel	0.23mm (32G) x 4mm	100
<b>IV505325</b>	Insulin Pen Needle 5 Bevel	0.23mm (32G) x 5mm	100
<b>IV505326</b>	Insulin Pen Needle 5 Bevel	0.23mm (32G) x 6mm	100