

Dual Safety Insulin Pen Needles

Our Dual Safety Insulin Pen Needles provide safety and comfort through their advanced technology.

The user is protected from accidental needlestick by the automatic safety shields, which activate at both ends, after injection.

The broad shield diameter offers greater comfort by reducing pressure on the skin and facilitating penetration.

The mdevices Dual Safety Insulin Pen Needles fit a broad variety of pens including; Sanofi Aventis, Lilly, Novo Nordisk, Owen Mumford and Berlin Chemie.

A choice of needle gauges (30G and 31G) and needle lengths (4mm, 5mm, 6mm and 8mm) are available and are coloured coded for easy identification.

- ⊗ Dual automatic shield lock for safety
- ⊗ Less resistance during needle penetration
- ⊗ Lubricated needle for a smooth injection
- ⊗ Several needle gauges and sizes available
- ⊗ Compatible with a broad variety of insulin pens
- ⊗ Packaged in a sealed sterility barrier



Automatic Safety Shield 1



Automatic Safety Shield 2

Dual Safety Insulin Pen Needles

Dual Safety Insulin Pen Needles				STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY
IV512305	Dual Safety Insulin Pen Needle	0.3mm (30G) x 5mm	Purple	100
IV512306	Dual Safety Insulin Pen Needle	0.3mm (30G) x 6mm	Orange	100
IV512308	Dual Safety Insulin Pen Needle	0.3mm (30G) x 8mm	Blue	100
IV512314	Dual Safety Insulin Pen Needle	0.25mm (31G) x 4mm	Green	100
IV512315	Dual Safety Insulin Pen Needle	0.25mm (31G) x 5mm	Purple	100
IV512316	Dual Safety Insulin Pen Needle	0.25mm (31G) x 6mm	Orange	100