


## Auto-retractable Safety Insulin Syringe

m|devices auto-retractable Safety Insulin Syringes provide enhanced safety and ease of use. The needle automatically retracts into the syringe barrel immediately following injection, significantly reducing the risk of accidental needle stick injuries.

One-handed activation is achieved through a simple, continuous press of plunger to its end.

- ⊗ Automatic needle retraction for reduced needle stick risk
- ⊗ One handed safety activation
- ⊗ Bold scale markings
- ⊗ Orange cap for easy identification
- ⊗ Pre-attached needle
- ⊗ Clear barrel
- ⊗ Peel open pouch



Safety Insulin Syringe - Auto Retractable			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
IV540298	Safety Insulin Syringes - Auto-retractable, 29G x 8mm	0.5mL	100
IV541298	Safety Insulin Syringes - Auto-retractable, 29G x 8mm	1mL	100
IV541308	Safety Insulin Syringes - Auto-retractable, 30G x 8mm	1mL	100