

Cervical Dilating Balloon Procedure Pack

Efficiency

Our 'all-in-one' improves performance by eliminating collection times.

Environmental Impact

Significantly reduces individual packaging waste

Fit for Purpose

All the necessary items are in one pack, specifically designed to support the insertion of m|devices Cervical Dilating Balloon

- ⊗ Enables faster more efficient setup
- ⊗ Eliminates time spent collecting individual items
- ⊗ Assists with maintaining infection control principles
- ⊗ Significantly reduces individual packaging waste
- ⊗ Specifically designed "fit" for purpose procedure pack
- ⊗ Inclusion of the m|devices Cervical Dilating Balloon 100% Silicone, closed distal end with high contrast position markings printed on the shaft
- ⊗ Includes batch tracking sticker
- ⊗ Latex free



Cervical Dilating Balloon Procedure Pack

STERILE

RE-ORDER	DESCRIPTION	CARTON QTY
OB010002	Cervical Dilating Balloon Procedure Pack	30

m|devices[®]

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321

Cervical Dilating Balloon Procedure Pack

Pack Items

QUANTITY	DESCRIPTION	REF.	QUANTITY	DESCRIPTION	REF.
1	Sterile Field 120cx120cm	A	1	m devices Cervical Dilating Balloon	G
1	Blue Rectangle Tray	B	1	XRD Gauze Swab 10x10cm 8Ply (5)	H
2	Lubricating Gel Sterile 27g	C	1	m devices Securement Device	I
4	Luer Lock Syringe 20mL	D	1	Plastic Forceps Clear	J
2	Gallipots 120mL	E	1	Tracking Label	K
1	SS Sponge Forceps 25cm	F			

