

## Cervical Dilating Balloon Procedure Pack

### Efficiency

Our 'all-in-one' improves performance by eliminating collection times.

### Environmental Impact

Significantly reduces individual packaging waste

### Environmental Impact

Significantly reduces individual packaging waste

- ⊗ Enables faster more efficient setup
- ⊗ Eliminates time spent collecting individual items
- ⊗ Assists with maintaining infection control principles
- ⊗ Significantly reduces individual packaging waste
- ⊗ Specifically designed "fit" for purpose procedure pack
- ⊗ Inclusion of the m|devices Cervical Dilating Balloon 100% Silicone, closed distal end with high contrast position markings printed on the shaft
- ⊗ Includes batch tracking sticker
- ⊗ Latex free



Cervical Dilating Balloon Procedure Pack		STERILE
RE-ORDER	DESCRIPTION	CARTON QTY
OB010002	Cervical Dilating Balloon Procedure Pack	30

**m|devices®**

**m|devices**  
27 Llewellyn Avenue  
Villawood NSW 2163  
Australia

T +61 (0)2 8718 2800  
E [info@mdevices.com.au](mailto:info@mdevices.com.au)  
W [mdevices.com.au](http://mdevices.com.au)

Cervical Dilating Balloon Procedure Pack

Pack Items					
QUANTITY	DESCRIPTION	REF.	QUANTITY	DESCRIPTION	REF.
1	Sterile Field 120cx120cm	A	1	m devices Cervical Dilating Balloon	G
1	Blue Rectangle Tray	B	1	XRD Gauze Swab 10x10cm 8Ply (5)	H
2	Lubricating Gel Sterile 27g	C	1	m devices Securement Device	I
4	Luer Lock Syringe 20mL	D	1	Plastic Forceps Clear	J
2	Gallipots 120mL	E	1	Tracking Label	K
1	SS Sponge Forceps 25cm	F			

