


Cervical Dilating Balloon Introducer

The Cervical Dilating Balloon Introducer has been designed in collaboration with clinicians to assist with the insertion of the m|devices Cervical Dilating Balloon.

Its malleable nylon material enables a simple and efficient insertion and the J-loop handle allows the placement of the Cervical Dilating Balloon with minimal movement.

The Introducer can be removed by simply reverse-bending the J-loop handle.

- ⊗ Bespoke design to aid the insertion of the Cervical Dilating Balloon
- ⊗ Malleable nylon material to support insertion
- ⊗ J-loop handle enables the introducer to be secured in place during use
- ⊗ Easily manipulated into any required angle or shape
- ⊗ Designed for easy insertion and removal



Used with the Cervical Dilating Balloon (sold separately)

Cervical Dilating Balloon Introducer

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
OB020000	Cervical Dilating Balloon Introducer	43.5cm	20

m|devices[®]

m|devices
27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800
E info@mdevices.com.au
W mdevices.com.au