



2-Way Foley Catheters


Approved for both urethral and suprapubic* catheterisation.

They feature a symmetrical balloon design to support reliable positioning and the unisex range comes with high contrast incremental markings to aid catheter placement.

- ⊗ Intended for Urethral and SPC use
- ⊗ 100% Silicone
- ⊗ High contrast incremental markings from 6cm to 30cm (on Unisex range)
- ⊗ Long indwelling time of up to 8 weeks**
- ⊗ Smooth rounded tip to aid a gentle insertion
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Three lengths available: Paediatric 33cm, Female 23cm and Unisex 45cm
- ⊗ Colour coded for easy size identification



Standard Tip



Incremental markings on Unisex range to support placement

2-Way Foley Catheters

STERILE

RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR010006	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 1.5mL Balloon	6Fr	Pink	10
UR010000	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon	8Fr	Light Blue	10
UR010001	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon	10Fr	Black	10
UR010012	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	12Fr	White	10
UR010014	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	14Fr	Green	10
UR010016	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	16Fr	Orange	10

Table continued on next page

m|devices®

m|devices
27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800
E info@mdevices.com.au
W mdevices.com.au

2-Way Foley Catheters

2-Way Foley Catheters				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR011000	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	12Fr	White	10
UR011001	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	14Fr	Green	10
UR011002	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	16Fr	Orange	10
UR011003	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	18Fr	Red	10
UR011004	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	20Fr	Yellow	10
UR011005	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	22Fr	Purple	10
UR011006	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	24Fr	Blue	10
UR011007	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	26Fr	Pink	10
UR011016	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	16Fr	Orange	10
UR011017	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	18Fr	Red	10
UR011018	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	20Fr	Yellow	10
UR011019	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	22Fr	Purple	10
UR011020	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	24Fr	Blue	10
UR011021	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	26Fr	Pink	10

* SPC use excludes paediatric range.

** As per laboratory-based testing results and validation against requirements ISO 20696.