


2-Way Foley Catheters

Approved for both urethral and suprapubic catheterisation, our 2-Way Foley is made from 100% medical grade translucent and flexible silicone. A symmetrical balloon design supports reliable positioning and minimises the potential for ridging.



- ⊗ 100% Silicone
- ⊗ Available in male, female and paediatric lengths
- ⊗ Intended for Urethral and SPC use
- ⊗ Long indwelling time of up to 8 weeks*
- ⊗ Colour coded for easy size identification
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Radiopaque

| 2-Way Foley Catheters | | | | | STERILE |
|-----------------------|--------------------------------------------------------------------------|-------|------------|------------|---------|
| RE-ORDER | DESCRIPTION | GAUGE | COLOUR | CARTON QTY | |
| UR010006 | 2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 1.5mL Balloon | 6Fr | Pink | 10 | |
| UR010000 | 2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon | 8Fr | Light Blue | 10 | |
| UR010001 | 2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon | 10Fr | Black | 10 | |
| UR010012 | 2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon | 12Fr | White | 10 | |
| UR010014 | 2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon | 14Fr | Green | 10 | |
| UR010016 | 2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon | 16Fr | Orange | 10 | |
| UR011000 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 12Fr | White | 10 | |
| UR011001 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 14Fr | Green | 10 | |
| UR011002 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 16Fr | Orange | 10 | |

Table continued on next page

2-Way Foley Catheters

| 2-Way Foley Catheters | | | | STERILE | |
|-----------------------|------------------------------------------------------------|-------|--------|------------|--|
| RE-ORDER | DESCRIPTION | GAUGE | COLOUR | CARTON QTY | |
| UR011003 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 18Fr | Red | 10 | |
| UR011004 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 20Fr | Yellow | 10 | |
| UR011005 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 22Fr | Purple | 10 | |
| UR011006 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 24Fr | Blue | 10 | |
| UR011007 | 2-Way Foley Catheter, Standard Tip, 40cm with 10mL Balloon | 26Fr | Pink | 10 | |
| UR011016 | 2-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 16Fr | Orange | 10 | |
| UR011017 | 2-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 18Fr | Red | 10 | |
| UR011018 | 2-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 20Fr | Yellow | 10 | |
| UR011019 | 2-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 22Fr | Purple | 10 | |
| UR011020 | 2-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 24Fr | Blue | 10 | |
| UR011021 | 2-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 26Fr | Pink | 10 | |