


3-Way Silicone Foley Catheters

The 3-Way Foley is made from 100% medical grade translucent and flexible silicone. The symmetrical balloon design supports reliable positioning in the bladder and minimises the potential for ridging.

- ⊗ Extended 43cm length
- ⊗ Long indwelling time of up to 29 days
- ⊗ Colour-coded to support easy size identification
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Radiopaque



3-Way Silicone Foley Catheters

STERILE

RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR013000	Silicone 3-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon	16Fr	Orange	10
UR013001	Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 30mL Balloon	18Fr	Red	10
UR013002	Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 50mL Balloon	20Fr	Yellow	10
UR013003	Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 50mL Balloon	22Fr	Purple	10
UR013004	Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 60mL Balloon	24Fr	Blue	10

m|devices

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321