


3-Way Silicone Foley Catheters

The 3-Way Foley is made from 100% medical grade translucent and flexible silicone. The symmetrical balloon design supports reliable positioning in the bladder and minimises the potential for ridging.

- ⊗ Extended 43cm length
- ⊗ Long indwelling time of up to 30 days
- ⊗ Colour-coded to support easy size identification
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Radiopaque



3-Way Silicone Foley Catheters

STERILE

| RE-ORDER | DESCRIPTION | GAUGE | COLOUR | CARTON QTY |
|----------|---|-------|--------|------------|
| UR013000 | Silicone 3-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon | 16Fr | Orange | 10 |
| UR013001 | Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 30mL Balloon | 18Fr | Red | 10 |
| UR013002 | Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 50mL Balloon | 20Fr | Yellow | 10 |
| UR013003 | Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 50mL Balloon | 22Fr | Purple | 10 |
| UR013004 | Silicone 3-Way Foley Catheter, Standard Tip, 43cm with 60mL Balloon | 24Fr | Blue | 10 |

m|devices

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321