


Irrigation Sets - Cystoscopy

Our Cystoscopy Irrigation Sets are made of high quality flexible tubing with a rigid roller regulator to accurately adjust flow rates according to clinical need during procedures.

- ⌵ High flow rates
- ⌵ Rigid roller regulator for controlled flow
- ⌵ Soft pillow drip chamber
- ⌵ Single spike
- ⌵ Stepped connector with silicone tube
- ⌵ Available with male luer lock connector



Irrigation Sets - Cystoscopy			STERILE
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
UR020000	Cystoscopy Irrigation Set	260cm	50
UR020002	Cystoscopy Irrigation Set with MLL	285cm	50