


Irrigation Set - Cystoscopy

This Cystoscopy TUR Irrigation Set is uniquely designed for high flow irrigation post transurethral resection of the prostate or bladder to support the removal of debris, clots and tissue. The flexible soft tubing and drip chamber, double spike configuration and stepped connector with silicone tube supports the full range of clinical needs.

- ⊗ High flow rates
- ⊗ Double spike configuration
- ⊗ Soft chamber without filter
- ⊗ Two click clamps
- ⊗ Stepped connector with silicone tube



Irrigation Set - Cystoscopy			STERILE
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
UR020001	Cystoscopy Irrigation Set	285cm	50