


## Hydrophilic Coated Nelaton Catheters

Our Hydrophilic catheters are made from non-toxic DEHP free PVC. Their smooth water soluble polymer coating combined with a heat polished distal tip supports a smooth insertion without the need of lubricants.

- ⊗ Available in male, female and paediatric lengths
- ⊗ Heat polished distal end aids insertion
- ⊗ Hydrophilic coating supports lubricant free insertion
- ⊗ Colour coded to support easy size identification
- ⊗ DEHP free



Hydrophilic Coated Nelaton Catheter				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR043001	Hydrophilic Coated Nelaton Catheter, 30cm (Paediatric)	8Fr	Light Blue	25
UR043002	Hydrophilic Coated Nelaton Catheter, 30cm (Paediatric)	10Fr	Black	25
UR044000	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	8Fr	Light Blue	25
UR044001	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	10Fr	Black	25
UR044002	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	12Fr	White	25

Table continued on next page

## Hydrophilic Coated Nelaton Catheters

Hydrophilic Coated Nelaton Catheter				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
<b>UR044004</b>	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	16Fr	Orange	25
<b>UR044005</b>	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	18Fr	Red	25
<b>UR045000</b>	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	8Fr	Light Blue	25
<b>UR045001</b>	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	10Fr	Black	25
<b>UR045002</b>	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	12Fr	White	25
<b>UR045003</b>	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	14Fr	Green	25
<b>UR045004</b>	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	16Fr	Orange	25
<b>UR045005</b>	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	18Fr	Red	25