


Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip

Our Hydrophilic Nelaton Catheter with Water Sachet- Coude Tip, is an 'all-in-one' product that is convenient and easy to use.

The angled tip makes it easy to navigate the urethra and supports a difficult insertion.

- ⊗ Supports a difficult insertion
- ⊗ No-touch sleeve
- ⊗ Easy peel packaging with wall attachment
- ⊗ Hydrophilic coating supports lubricant free insertion
- ⊗ Heat polished distal end aids insertion



Hydrophilic Coated Nelaton Catheter with Water Sachet - Coude Tip

STERILE

RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR049000	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	8Fr	Light Blue	25
UR049001	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	10Fr	Black	25
UR049002	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	12Fr	White	25
UR049003	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	14Fr	Green	25
UR049004	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	16Fr	Orange	25
UR049005	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	18Fr	Red	25