


## Intermittent Closed Catheterisation System

Our Closed Catheterisation System is sterile and pre-lubricated with a hydrophilic gel, so no additional preparation is needed. The system works like a traditional Nelaton Catheter and provides a simple catheterisation that reduces the risk of urinary tract infections.

- ⊗ Integrated 1500mL bag
- ⊗ Sterile closed system reduces the chance of infection
- ⊗ Pre-lubricated for simple insertion
- ⊗ DEHP, BPA and latex free



Intermittent Closed Catheterisation System			STERILE
RE-ORDER	DESCRIPTION	GAUGE	CARTON QTY
UR050008	Standard Intermittent Catheter with Gel and 1500mL Bag	8Fr	25
UR050010	Standard Intermittent Catheter with Gel and 1500mL Bag	10Fr	25
UR050012	Standard Intermittent Catheter with Gel and 1500mL Bag	12Fr	25
UR050014	Standard Intermittent Catheter with Gel and 1500mL Bag	14Fr	25
UR050016	Standard Intermittent Catheter with Gel and 1500mL Bag	16Fr	25