

Urine Bags

The m|devices 2500mL hospital Urine Bag is designed to meet clinicians' needs with an integrated reinforced double hanger, rope hanger and a bed sheet clamp to secure the bag. The 120cm tube length has an integrated needleless sampling port. The distal tube end, beneath the bag, has a T-Tap and it also features an integrated anti-reflux valve.


- 2500mL Bag
- Multiple securement options; integrated hanger, rope hanger and bedsheet clamp
- 120cm tubing with needleless sampling port
- T-Tap with protective housing to minimise infection
- Anti-reflux valve to prevent backflow to reduce infection
- Clear incremental measurement markings



Bedsheet clamp



Needleless sampling port



T-Tap

Urine Bag

STERILE

RE-ORDER	DESCRIPTION	BAG CAPACITY	TUBING LENGTH	CARTON QTY
UR062500	Hospital Urine Bag with NRV and Integrated Hanger	2500mL	120cm	40

m|devices[®]

m|devices
27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800
E info@mdevices.com.au
W mdevices.com.au