

## Urine Bag Hanger

Our easy-to-use reinforced double Urine Bag Hanger features prongs to attach the bag and also includes a clasp to secure the urine bag tubing in place.

- Suitable for all major bag types
- Reinforced for effective securement
- Easy to use prongs to attach the urine bag
- Pre-formed securement clasp to secure the urine bag tubing in place



### Urine Bag Hanger

NON-STERILE

RE-ORDER DESCRIPTION

CARTON QTY

UR063003NS Urine Bag Hanger

50

**m|devices**<sup>®</sup>

**m|devices**  
27 Llewellyn Avenue  
Villawood NSW 2163  
Australia

T +61 (0)2 8718 2800  
E [info@mdevices.com.au](mailto:info@mdevices.com.au)  
W [mdevices.com.au](http://mdevices.com.au)