

## Leg Bags - T-Tap

Our Leg Bags are designed to be unobtrusive offering discreet storage under clothing. They have a T-Tap and non-return valve to support a sterile fluid pathway. The soft non-woven backing improves user comfort and they are available with either an adjustable elastic fabric strap or silicone lined strap.



- ⊗ Discreet to wear
- ⊗ Soft non-woven backing is comfortable next to the skin
- ⊗ Adjustable elastic fabric strap or silicone lined strap
- ⊗ Multiple tube lengths; 10cm, 30cm and 50cm
- ⊗ Various volume capacities; 350mL, 500mL and 750mL
- ⊗ Bonded step connector
- ⊗ The T-Tap design prevents accidental drainage sometimes experienced with other tap styles



| Leg Bags  |                                                                |          | STERILE    |
|-----------|----------------------------------------------------------------|----------|------------|
| RE-ORDER  | DESCRIPTION                                                    | CAPACITY | CARTON QTY |
| UR135010* | LEG BAG - T-Tap NRV 10cm Short Tube with Silicone Lined Straps | 350mL    | 80         |
| UR135030* | LEG BAG - T-Tap NRV 30cm Long Tube with Silicone Lined Straps  | 350mL    | 80         |
| UR135050  | LEG BAG - T-Tap NRV 50cm Long Tube with Silicone Lined Straps  | 350mL    | 80         |

Table continued on next page

**m|devices**

m|devices  
27 Llewellyn Avenue  
Villawood NSW 2163  
Australia

T +61 (0)2 8718 2800  
E info@mdevices.com.au  
W mdevices.com.au

## Leg Bags - T-Tap

| Leg Bags         |                                                                |          | STERILE    |
|------------------|----------------------------------------------------------------|----------|------------|
| RE-ORDER         | DESCRIPTION                                                    | CAPACITY | CARTON QTY |
| <b>UR150010*</b> | LEG BAG - T-Tap NRV 10cm Short Tube with Silicone Lined Straps | 500mL    | 80         |
| <b>UR150030</b>  | LEG BAG - T-Tap NRV 30cm Long Tube with Silicone lined Straps  | 500mL    | 80         |
| <b>UR150050*</b> | LEG BAG - T-Tap NRV 50cm Long Tube with Silicone lined Straps  | 500mL    | 80         |
| <b>UR175010</b>  | LEG BAG - T-Tap NRV 10cm Short Tube with Silicone lined Straps | 750mL    | 80         |
| <b>UR175030</b>  | LEG BAG - T-Tap NRV 30cm Long Tube with Silicone lined Straps  | 750mL    | 80         |
| <b>UR175050</b>  | LEG BAG - T-Tap NRV 50cm Long Tube with Silicone Lined Straps  | 750mL    | 80         |
| <b>UR235010*</b> | LEG BAG - T-Tap NRV 10cm Short Tube with Elastic Straps        | 350mL    | 100        |
| <b>UR235030*</b> | LEG BAG - T-Tap NRV 30cm Long Tube with Elastic Straps         | 350mL    | 100        |
| <b>UR235050*</b> | LEG BAG - T-Tap NRV 50cm Long Tube with Elastic Straps         | 350mL    | 100        |
| <b>UR250010</b>  | LEG BAG - T-Tap NRV 10cm Short Tube with Elastic Straps        | 500mL    | 100        |
| <b>UR250030</b>  | LEG BAG - T-Tap NRV 30cm Long Tube with Silicone Lined Straps  | 500mL    | 100        |
| <b>UR250050*</b> | LEG BAG - T-Tap NRV 50cm Long Tube with Elastic Straps         | 500mL    | 100        |
| <b>UR275010</b>  | LEG BAG - T-Tap NRV 10cm Short Tube with Elastic Straps        | 750mL    | 100        |
| <b>UR275030*</b> | LEG BAG - T-Tap NRV 30cm Long Tube with Elastic Straps         | 750mL    | 100        |
| <b>UR275050*</b> | LEG BAG - T-Tap NRV 50cm Long Tube with Elastic Straps         | 750mL    | 100        |

\* This product is currently being sold without the tapered shape. All stock will transition to the tapered shape upon depletion of existing stock.