


Leg Bags - Lever Tap

Our Leg Bags are designed to be unobtrusive offering discreet storage under clothing. They have a Lever Tap design to support users with reduced manual dexterity. To improve comfort, they feature a soft non-woven backing, silicone lined strap and fabric fastener.



- ⊗ Discreet to wear
- ⊗ Soft non-woven backing is comfortable next to the skin
- ⊗ Silicone lined straps and fabric fastener
- ⊗ Multiple tube lengths; 10cm, 30cm and 50cm
- ⊗ Various volume capacities; 350mL, 500mL and 750mL
- ⊗ Bonded step connector
- ⊗ Supports users who have reduced manual dexterity as the Lever Tap can be easily be opened and closed



Silicone lined strap



Lever Tap

Leg Bags			STERILE
RE-ORDER	DESCRIPTION	CAPACITY	CARTON QTY
UR335010	LEG BAG - Lever Tap NRV 10cm Short Tube with Silicone lined Straps	350mL	80
UR335030	LEG BAG - Lever Tap NRV 30cm Long Tube with Silicone lined Straps	350mL	80
UR335050	LEG BAG - Lever Tap NRV 50cm Long Tube with Silicone lined Straps	350mL	80

Table continued on next page

m|devices

m|devices
27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800
E info@mdevices.com.au
W mdevices.com.au

Leg Bags - Lever Tap

Leg Bags			STERILE
RE-ORDER	DESCRIPTION	CAPACITY	CARTON QTY
UR350010	LEG BAG - T-Tap NRV 10cm Short Tube with Silicone Lined Straps	500mL	80
UR350030	LEG BAG - T-Tap NRV 30cm Long Tube with Silicone Lined Straps	500mL	80
UR350050	LEG BAG - T-Tap NRV 50cm Long Tube with Silicone Lined Straps	500mL	80
UR375010	LEG BAG - T-Tap NRV 10cm Short Tube with Silicone lined Straps	750mL	80
UR375030	LEG BAG - T-Tap NRV 30cm Long Tube with Silicone lined Straps	750mL	80
UR375050	LEG BAG - T-Tap NRV 50cm Long Tube with Silicone Lined Straps	750mL	80