

Yeescope™

Disposable Laryngoscope

PRODUCT DATA SHEET



DESCRIPTION

m|devices Yeescope™ is a one-piece disposable laryngoscope with incorporated handle and blade. It has a light weight; patented angle design has a bright white light that allows for powerful illumination. Single use design avoids cross contamination and incompatibility between blades and handles. Simple function testing of the light within the packaging makes this product perfect for resuscitation trolleys which require daily checks.

FEATURES AND BENEFITS

- Available in 2 sizes: MAC3 and MAC4
- Light weight design
- Patented design
- Easy illumination activation
- Prevents cross contamination
- Eliminates malfunction due to handle and blade compatibility

PRODUCT DATA

- Battery (regular): AA
- MAC3 proven suitable: 3mths to adults
- Sterile EO
- Latex free
- DEHP: Yes
- Material: ABS
- Spanner available: AN010002
- Shelf Life: 3 Years
- ARTG Number (regular): 210822

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
AN010000	Yeescope - Fully Disposable Laryngoscope	MAC3	10
AN010001	Yeescope - Fully Disposable Laryngoscope	MAC4	10