

Laryngeal Mask Airways

PRODUCT DATA SHEET



DESCRIPTION

m|devices Laryngeal Mask Airways are made from 100% medical grade silicone. A smooth one-piece design with guide markings on the handle for easy insertion and positioning. The translucent tubing allows clear visibility and the soft seal cuff which is permeable to nitrous oxide, reduces pressure rise and minimises potential trauma.

FEATURES AND BENEFITS

- 100% Medical grade silicone
- Smooth one-piece design
- Guide markings
- Colour coded for easy size selection
- Translucent tubing
- Soft seal cuff
- Peel open pouch
- Single patient use only

PRODUCT DATA

- Sterile EO
- Latex free
- DEHP free
- Material: Silicone
- Shelf Life: 3 Years
- ARTG Number: 210988

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
AN030000	Laryngeal Mask Airway	1	10
AN030001	Laryngeal Mask Airway	1.5	10
AN030002	Laryngeal Mask Airway	2	10
AN030003	Laryngeal Mask Airway	2.5	10
AN030004	Laryngeal Mask Airway	3	10
AN030005	Laryngeal Mask Airway	4	10
AN030006	Laryngeal Mask Airway	5	10