

# Suction Catheters

## Open Tip Funnel Connector

### DESCRIPTION

The transparent m|devices Suction Catheter with Funnel Connectors are designed for easy passage through endotracheal or tracheostomy tubing. Available in paediatric through to adult sizes, with two different lengths and a variety of diameters. The Suction Catheters are made from clear and soft non-toxic medical grade PVC. Smooth open distal tip with an elliptical eye for easy and effective suctioning. Additional inner sleeve with tear off perforations to assist the clinician with maintaining aseptic technique whilst inserting the Suction catheter.

### FEATURES & BENEFITS

- Funnel connector allows for a secure connection into suction tubing
- Wide range of sizing options available
- Colour coded to assist with easy size identification
- Inner sleeve packaging with perforations to assist with maintaining aseptic technique during insertion
- Easy to open peel pouch
- Made from transparent PVC to allow for visualisation when suctioning


### PRODUCT DATA

- Length: 380mm for 6fr and 500mm for sizes 8Fr-16Fr
- Size range: 6fr - 16fr
- Not made with DEHP
- Not made with Natural Rubber Latex
- Sterile EO
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228232

### STANDARDS

- ISO 13485:2003 Medical devices-Quality Management Systems - Requirements for Regulatory Purposes
- ENISO 14155-2:2003 Clinical investigation for medical devices for human subjects - part2: Clinical investigation plans
- EN550:1994 Sterilisation of medical devices - Validation and routine control of ethylene oxide sterilisation

### PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Funnel Connector
2	Open Tip

# Suction Catheters

## Open Tip Funnel Connector

Suction Catheters Open Tip (Funnel Connector)					
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
<b>AN041000</b>	Suction Catheter, Open Tip, Funnel Connector	380mm	6Fr	Light Green	50
<b>AN041001</b>	Suction Catheter, Open Tip, Funnel Connector	500mm	8Fr	Light Blue	50
<b>AN041002</b>	Suction Catheter, Open Tip, Funnel Connector	500mm	10Fr	Black	50
<b>AN041003</b>	Suction Catheter, Open Tip, Funnel Connector	500mm	12Fr	White	50
<b>AN041004</b>	Suction Catheter, Open Tip, Funnel Connector	500mm	14Fr	Green	50
<b>AN041005</b>	Suction Catheter, Open Tip, Funnel Connector	500mm	16Fr	Orange	50