

# Suction Catheters

Open Tip  
Y Type Control Vent

## DESCRIPTION

The transparent m|devices Suction Catheters, Open Tip with Y Type Control Vent, are designed for easy passage through endotracheal or tracheostomy tubing. Available in a variety of diameters for both paediatric and adult use with lengths of 100mm, 380mm and 560mm. The catheters are made from clear non-toxic medical grade PVC. Smooth open distal end with elliptical eye for easy and effective suctioning. Additional inner sleeve with tear off perforations to assist the clinician with maintaining aseptic technique whilst inserting the Suction Catheter.

Medium style packaging is available to allow for easy storage in space limited areas.

## FEATURES & BENEFITS

- Stepped connector for easy connection to suction tubing
- Wide range of size options
- Incremental measurement markings
- Inner sleeve with perforations along to assist maintaining aseptic technique
- Colour coded for ease of size selection
- Peel open pouch
- Made from translucent PVC


## PRODUCT DATA

- Length: 380mm, 560mm
- Size range: 5Fr-18Fr
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Not made with Natural Rubber Latex
- Sterile EO
- Material: PVC
- Shelf Life: 5 Years
- ARTG Number: 228232

## STANDARDS

- ISO 13485:2003 Medical devices-Quality Management Systems - Requirements for Regulatory Purposes
- ENISO 14155-1:2003 Clinical investigation for medical devices for human subjects - part1: General Requirements
- ENISO 14155-2:2003 Clinical investigation for medical devices for human subjects - part2: Clinical investigation plans
- EN550:1994 Sterilisation of medical devices - Validation and routine control of ethylene oxide sterilisation

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Y Type Control Vent
2	Colour Coded
3	Open Tip

# Suction Catheters

## Open Tip Y Type Control Vent

Suction Catheters Open Tip (Y Type Control Vent)					
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
<b>AN045002</b>	Suction Catheter, Open Tip, Y Type Control Vent	100mm	8Fr	Light Blue	50
<b>AN045003</b>	Suction Catheter, Open Tip, Y Type Control Vent	100mm	10Fr	Black	50
<b>AN045004</b>	Suction Catheter, Open Tip, Y Type Control Vent	100mm	12Fr	White	50
<b>AN042000*</b>	Suction Catheter, Open Tip, Y Type Control Vent	380mm	5Fr	Grey	50
<b>AN042001*</b>	Suction Catheter, Open Tip, Y Type Control Vent	380mm	6Fr	Light Green	50
<b>AN042008*</b>	Suction Catheter, Open Tip, Y Type Control Vent	380mm	7Fr	Ivory	50
<b>AN042002</b>	Suction Catheter, Open Tip, Y Type Control Vent	560mm	8Fr	Light Blue	50
<b>AN042003</b>	Suction Catheter, Open Tip, Y Type Control Vent	560mm	10Fr	Black	50
<b>AN042004</b>	Suction Catheter, Open Tip, Y Type Control Vent	560mm	12Fr	White	50
<b>AN042005</b>	Suction Catheter, Open Tip, Y Type Control Vent	560mm	14Fr	Green	50
<b>AN042006</b>	Suction Catheter, Open Tip, Y Type Control Vent	560mm	16Fr	Orange	50
<b>AN042007</b>	Suction Catheter, Open Tip, Y Type Control Vent	560mm	18Fr	Red	50
<b>AN046000*</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	380mm	5Fr	Grey	50
<b>AN046001*</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	380mm	6Fr	Light Green	50
<b>AN046002</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	8Fr	Light Blue	50
<b>AN046003</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	10Fr	Black	50
<b>AN046004</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	12Fr	White	50
<b>AN046005</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	14Fr	Green	50
<b>AN046006</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	16Fr	Orange	50
<b>AN046007</b>	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	18Fr	Red	50

\* 5-7Fr Suction Catheters do not feature the polished eyelets available on larger sizes.