

# Suction Tubing - Non-Sterile

## DESCRIPTION

The m|devices transparent non-sterile suction tubing is available in three versions and lengths; soft, flexible and ribbed. Standard hexangular connectors are attached to both ends, with some codes from the soft range that have options for only one hexangular connector. The tubing is made from non-toxic, medical grade PVC.

## FEATURES & BENEFITS

- Special soft, flexible or ribbed tubing
- Wide range of length options
- Tubing diameter:  
Flexible ID6mm OD9mm  
Rib ID6mm OD9mm
- Non-sterile
- With hexangular connectors AN054002 & AN054006  
Rib ID6.8mm OD10mm  
Soft ID6.8mm OD10mm
- Yellow coloured connectors to indicate Non-Sterile


## PRODUCT DATA

- Not made with DEHP
- Not Made with Natural Rubber Latex
- Medical Grade PVC
- ARTG Number: 211736
- Shelf Life: 5 Years
- Independant testing exceeds minimum requirements for suction tubing degree of collapse as per AS 2120.3 Clauses 5.2.2

## STANDARDS

- ISO 14971:2007 Medical devices-Application of risk management to medical devices
- ISO10993-1:2003 Biological evaluation of the medical devices-Part 1: Guidance on selection of tests
- ISO11135-1:2007 Sterilization of health care products-Ethylene oxide-Part1:Requirements for development, validation and routine control of a sterilisation process for medical devices

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tubing
2	Connector

# Suction Tubing - Non-Sterile

Suction Tubing - Non-Sterile					
RE-ORDER	DESCRIPTION	TUBING LENGTH (m)	TUBING TYPE	CARTON/QTY	
<b>AN050000YNS</b>	Suction Tubing - ID6mm OD9mm	2m	Flexible	80	
<b>AN050001YNS</b>	Suction Tubing - ID6mm OD9mm	3	Flexible	65	
<b>AN050002YNS</b>	Suction Tubing - ID6mm OD9mm	4.5m	Flexible	45	
<b>AN050003YNS</b>	Suction Tubing - ID6mm OD9mm	6m	Flexible	35	
<b>AN050004YNS</b>	Suction Tubing - ID6mm OD9mm	8m	Flexible	30	
<b>AN051000YNS</b>	Suction Tubing - ID6mm OD9mm	2m	Ribbed	80	
<b>AN051001YNS</b>	Suction Tubing - ID6mm OD9mm	3m	Ribbed	65	
<b>AN051002YNS</b>	Suction Tubing - ID6mm OD9mm	4.5m	Ribbed	45	
<b>AN051003YNS</b>	Suction Tubing - ID6mm OD9mm	6m	Ribbed	35	
<b>AN051004YNS</b>	Suction Tubing - 2 Connectors	8m	Ribbed	30	
<b>AN056001YNS</b>	Suction Tubing - Special Soft	3m	Special Soft	65	
<b>AN056006YNS</b>	Suction Tubing - Special Soft	4m	Special Soft	50	
<b>AN057001YNS</b>	Suction Tubing - Special Soft	3m	Special Soft	50	
<b>AN059000NS</b>	Suction Tubing - Special Soft with stepped connector	0.48m	Special Soft	150	