

Suction Tubing - Non-Sterile

DESCRIPTION

The m|devices transparent non-sterile suction tubing is available in three versions and lengths; soft, flexible and ribbed. Standard hexangular connectors are attached to both ends, with some codes from the soft range that have options for only one hexangular connector. The tubing is made from non-toxic, medical grade PVC.

FEATURES & BENEFITS

- Special soft, flexible or ribbed tubing
- Wide range of length options
- Tubing diameter:
Flexible ID6mm OD9mm
Rib ID6mm OD9mm
- Non-sterile
- With hexangular connectors AN054002 & AN054006
Rib ID6.8mm OD10mm
Soft ID6.8mm OD10mm
- Yellow coloured connectors to indicate Non-Sterile


PRODUCT DATA

- Not made with DEHP
- Not Made with Natural Rubber Latex
- Medical Grade PVC
- ARTG Number: 211736
- Shelf Life: 5 Years
- Independant testing exceeds minimum requirements for suction tubing degree of collapse as per AS 2120.3 Clauses 5.2.2

STANDARDS

- ISO 14971:2007 Medical devices-Application of risk management to medical devices
- ISO10993-1:2003 Biological evaluation of the medical devices-Part 1: Guidance on selection of tests
- ISO11135-1:2007 Sterilization of health care products-Ethylene oxide-Part1:Requirements for development, validation and routine control of a sterilisation process for medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tubing
2	Connector

Suction Tubing - Non-Sterile

Suction Tubing - Non-Sterile				
RE-ORDER	DESCRIPTION	TUBING LENGTH (m)	TUBING TYPE	CARTON/QTY
AN050000YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	2m	Flexible	80
AN050001YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	3	Flexible	65
AN050002YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	4.5m	Flexible	45
AN050003YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	6m	Flexible	35
AN050004YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	8m	Flexible	30
AN051000YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	2m	Ribbed	80
AN051001YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	3m	Ribbed	65
AN051002YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	4.5m	Ribbed	45
AN051003YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	6m	Ribbed	35
AN051004YNS	Suction Tubing - ID6mm OD9mm, 2 connectors	8m	Ribbed	30
AN056001YNS	Suction Tubing - ID6mm OD10mm, 1 connector	3m	Special Soft	65
AN056006YNS	Suction Tubing - ID6mm OD10mm, 1 connector	4m	Special Soft	50
AN057001YNS	Suction Tubing - ID6mm OD10mm, 2 connectors	3m	Special Soft	50
AN059000NS	Suction Tubing - Special Soft with stepped connector	0.48m	Special Soft	150